



# **A Systematic Review of Deep Learning Methods for ACL Tear Detection and Grading**

Laya Hamehkhani | Dr. Mahdi Hossieni  
Nov. 2025

# Table of content

- Introduction
- Segmentation-Based Approaches
- Detection & Localization Systems
- Multi-Class Grading & Severity Assessment
- Advanced Architecture Approaches
- Data-Efficient & Few-Shot Learning
- Multi-Plane & Multi-Task Models

# ACL Injuries - Clinical Burden

## The Clinical Challenge

- Epidemiology & Impact
  - **200,000+ ACL injuries** annually in the United States alone
  - Peak incidence: **16-39 years** (active population)
  - **Female** athletes: **2-8× higher risk** than **males** in comparable sports
- Diagnostic Complexity
  - Accuracy varies by reader experience: **86-94% (expert)** vs. **50-75% (general radiologist)**
  - Partial tears: Most challenging to diagnose, highest inter-reader variability
- Clinical Significance
  - Delays in diagnosis → Progressive joint instability → **Secondary injuries**
  - **Misclassification** of partial vs. complete tears → **Inappropriate treatment decisions**
  - Long-term: **50%** develop osteoarthritis within **10-20 years post-injury**

# Current Diagnostic Workflow & Bottlenecks

## MRI-Based Diagnosis: Standard of Care

- MRI Protocol for ACL Assessment
  - **Primary sequences:** Proton density (PD), T2-weighted, Fat-suppressed
  - **Key imaging planes:** Sagittal (primary), Coronal, Axial
  - **Target structures:** ACL morphology, signal intensity, fiber continuity
  - **Indirect signs:** Bone bruising, PCL angle, anterior tibial translation
- **Current Workflow Limitations**

Challenge	Impact	Quantification
<b>Reading time</b>	Radiologist workload	10-20 min per knee MRI
<b>Manual segmentation</b>	Research/surgical planning bottleneck	1-2 hours per case
<b>Inter-reader variability</b>	Inconsistent reporting	Cohen's $\kappa$ : 0.42-0.65 (moderate)
<b>Subtle tear detection</b>	Missed partial tears	15-30% miss rate for grade I-II
<b>Multi-abnormality assessment</b>	Cognitive overload	12+ potential pathologies to evaluate
<b>Experience dependency</b>	Quality variation across centers	20-40% accuracy gap

# Why Deep Learning for ACL Diagnosis?

The Convergence of Technology and Clinical Need

- Unique Advantages of DL in Musculoskeletal Imaging
  - 1. Pattern Recognition at Scale
    - Learns hierarchical features from raw pixel data
    - Captures subtle patterns invisible to human eye
    - Processes multi-planar, multi-sequence data simultaneously
  - 2. Consistency & Reproducibility
    - Eliminates inter-reader variability
    - Provides objective, quantifiable metrics
  - 3. Computational Efficiency
    - Real-time processing: <1 second vs. minutes
    - Scalable to high-volume clinical settings
    - Enables comprehensive multi-structure analysis
- From Research to Clinical Reality
  - FDA clearances emerging for orthopedic AI tools (2019-present)
  - Multiple commercial platforms in development
  - Integration with PACS/radiology workflows advancing

# MRI Fundamentals Quick Reference

## MRI Fundamentals Quick Reference

- MRI Sequences Mentioned in This Review

Sequence	Appearance	Primary Use in ACL Imaging
T1-weighted	Fat: bright, Fluid: dark	Anatomical detail, bone marrow
T2-weighted	Fluid: bright, Fat: intermediate	Edema, effusion detection
PD (Proton Density)	Balanced signal	ACL fiber visualization (primary)
Fat-Suppressed (FS)	Fat: dark, Fluid: very bright	Tear/edema sensitivity ↑

- Image Planes
  - **Sagittal:** Side-to-side slices (divides left/right)
  - **Coronal:** Front-to-back slices (divides front/back)
  - **Axial:** Top-to-bottom slices (divides superior/inferior)
- ACL MRI Appearance
  - **Normal:** Hypointense (dark), continuous fibers, 5-8mm thick, taut band
  - **Torn:** Discontinuous, hyperintense signal, irregular/absent fibers, wavy/lax

*Note: Deep learning models learn these patterns directly from data*

# Category 1: Segmentation-Based Approaches

## Precise Anatomical Structure Identification


- **Focus:** Automated pixel-level segmentation of knee anatomical structures including **ACL**, **PCL**, cartilage, meniscus, and surrounding tissues
- **Key Applications:**
  - **Pre-surgical planning and 3D visualization**
  - Reducing manual segmentation time (from 1-2 hours to <1 second)
  - **Multi-tissue analysis for comprehensive knee assessment**
  - Handling class imbalance challenges (small structures like ACL)
- **Clinical Value:**
  - Enables **precise measurement** of anatomical structures
  - Supports **holographic 3D modeling for surgical planning**
  - Foundation for downstream classification tasks

# HASA-ResUNet

## Hybrid Attention & Residual U-Net for Multi-Class Segmentation


- **Architecture:**

- Novel U-shaped encoder-decoder network
- Residual blocks replace standard convolutional layers
- Mitigate vanishing gradient problem



- **Key Innovations:**


- Addresses **pixel imbalance** (small structures like **ACL**)
- Handles **blurred boundaries** effectively
- Manages **shape diversity** in multi-class knee MRI segmentation



Visual segmentation results of the knee joint generated by the models.

# HASA-ResUNet

Hybrid Attention & Residual U-Net for Multi-Class Segmentation



- **HFEF Module:** Hierarchical feature fusion in skip connections → combines local + global features
- **ASA Module:** Atrous squeeze attention → multi-scale focus + long-range dependencies

# 3D DenseVNet

## Multi-Sequence 3D Semantic Segmentation

- **Architecture:**
  - **3D sampling pattern with encoder-decoder structure**
  - Trained on **multiple pulse sequences** as input channel: **T1 TSE, PD TSE, PD FS TSE, Angio GE**
- **Capabilities:**
  - Automatic segmentation of **13 distinct anatomical classes**:
    - PCL, ACL, cartilage, tendons, arteries, nerves, muscle
- **Key Finding:**
  - Best performance with all **four sequences** (T1PDFSAngio)
  - Confirms **multiple weightings enhance feature learning**
- **Clinical Applications:**
  - **Post & Pre-surgical visualization** (tibial drill holes, intra-articular grafts)
  - **3D rendering and holographic modeling**
  - Complete knee joint visualization with **transparent adipose tissue**

# 3D DenseVNet


## Multi-Sequence 3D Semantic Segmentation

(a) Post-surgical proton density fat suppressed TSE (FS TSE) sagittal plane image showing the torn ACL




(a)

(b) Predicted output by the neural network. The tibial drill hole and the intra-articular graft from the reconstructive surgery is visible.



(b)

(c) 3D render of the segmentation



(c)

(d) Holographic 3D model representation




(d)

# 3D DenseVNet

## Multi-Sequence 3D Semantic Segmentation



(a)



(b)



(c)



(d)



(e)




(f)



(g)



(h)




(i)

Results of test dataset in 3D

- (a), (b) Bone, ACL, PCL, meniscus and collateral ligaments
- (c) Bone, muscles, ligaments and veins
- (d) Complete segmentation with transparent adipose tissue
- (e) Cartilage
- (f) ACL
- (g) PCL
- (h), (i) Meniscus

# 3D DenseVNet

Multi-Sequence 3D Semantic Segmentation




The neural network architecture

# Modified 2D U\_Net FCNN


Ultra-Fast Automated Segmentation

- **Architecture:**
  - FCNN based on U-Net
- **Performance Breakthrough:**
  - Manual segmentation: **1-2 hours per subject**
  - Automated segmentation: **0.33 seconds per subject**
- **Impact:**
  - Eliminates major bottleneck (time-consuming process of segmentation) in clinical workflow
  - Enables real-time segmentation during MRI reading



# Modified 2D U\_Net FCNN

Ultra-Fast Automated Segmentation



Example modified 2D U-Net architecture (note: final model depth was tuned as hyperparameter).

# Conditional Generative Adversarial Networks (cGANs)


## Adversarial Learning for Multi-Tissue Segmentation

- **Architecture** (Based on **pix2pix** framework):
  - **Generator (G)**: Modified U-Net encoder-decoder
  - **Discriminator (D)**: PatchGAN (patch-based FCN)
- **Loss Function**:
  - cGAN loss + pixel-wise error loss
  - Ensures **output fidelity** and **reduces blurring**
- **Segmentation Targets**:
  - Femoral, tibial, and patellar bones
  - Cartilage surfaces
  - Cruciate ligaments
  - Selective muscles: medial vastus and gastrocnemius
- **Advantage**: Generates realistic, high-fidelity segmentation masks

Source Image




PatchGAN output with 1\*1 field size



# Conditional Generative Adversarial Networks (cGANs)

Adversarial Learning for Multi-Tissue Segmentation



Overview of cGANs architecture

# Category 2: Detection & Localization Systems

Pinpointing ACL Tears with Precision


- **Focus:** Identifying ACL **presence**, **location**, and **precise tear sites** within the ligament structure

- **Key Applications:**

- **Slice selection** (filtering ACL-containing images)
- **ROI isolation** (narrowing focus to intercondylar notch)
- **Rupture site classification** (femoral side, middle, tibial side)
- Automated **landmark detection**

- **Clinical Value:**



- Reduces radiologist workload by **filtering irrelevant images**
- **Provides surgical guidance for graft placement**
- Enables **precise documentation of injury location**



# Cascaded Detection with 3 CNNs

## Modular Three-Stage Pipeline

- **Three-Stage Process:**
  - Slice Detection (LeNet): Identifies ACL-containing images
  - Ligament Isolation (YOLO): Narrows to intercondylar notch ROI
  - Classification (DenseNet): Determines intact vs. torn
- **Design Philosophy:**
  - Highly modular cascaded structure
  - Emphasizes **precise localization before classification**
  - **Filters superfluous data** to enhance recognition accuracy
- **Outputs:**
  - Diagnosis result (Intact/Torn ACL)
  - Diagnosis probability map




Example of original image and its probability map of an ACL tear

# Cascaded Detection with 3 CNNs

Modular Three-Stage Pipeline

Illustration of the CNN architecture for the deep learning-based ACL tear detection system



# 3D CNN Scheme for Rupture Localization

Two-Step Segmentation and Landmark Detection

- **Step 1: Coarse Localization (Segmentation)**

- 3D U-Net architecture
- Categorizes knee into 4 components:
  - Femur, tibia, femoral footprints, tibial footprints


- **Step 2: Fine Localization & Classification**

- 3D Full CNN based on YOLOF (single-level detector)
- ResNet-101 backbone
- Heatmap regression network (adapted from 3D nnU-Net)

- **Classification Method:**

- Projects rupture point onto **line connecting footprint centers**
- Divides into **three equal parts**: femoral side, middle, tibial side


Also includes 2D CNN scheme with slice selection and landmark localization



An example of the model output

# 3D CNN Scheme for Rupture Localization

Two-Step Segmentation and Landmark Detection



# Category 3: Multi-Class Grading & Severity Assessment

Comprehensive WORMS-Based Staging (Whole-Organ Magnetic Resonance Imaging Score)

- **Focus:** Moving beyond **binary classification** to **multi-level severity grading** across multiple knee structures
- **Key Applications:**
  - WORMS-based severity staging (cartilage, bone, meniscus, ACL)
  - **Grade-specific treatment planning**
  - Standardized reporting aligned with clinical scales
  - Interreader agreement improvement
- **Clinical Value:**
  - Provides **nuanced diagnostic information**
  - Improves communication between radiologists
  - Supports **evidence-based treatment decisions**
  - **Confidence visualization reduces radiologist overload**

# Hierarchic 3D CNN Grading

## Multi-Class WORMS Severity Staging with AI Assistance


- **Scope:**
  - Full-knee 3D MRI processing
  - Hierarchic 3D Convolutional Classifiers
  - WORMS-based staging
- **Clinical Impact:**
  - AI boosts interreader agreement:
    - Cohen's  $\kappa$ : **0.42 (weak)  $\rightarrow$  0.61 (moderate)**
  - GUI displays confidence levels (**color + transparency**)

Colors indicate the lesion class:

red: severe | orange: moderate (or mild) | green: no lesion.

**Transparency:** the probability output by the model for that class (stronger colors  $\rightarrow$  more confidence)

The **ACL**, when deemed reconstructed, is set to a differentiated color with **no transparency level** and is indicated as **being postsurgical**.




# Hierarchic 3D CNN Grading

Multi-Class WORMS Severity Staging with AI Assistance

## Pipeline (4 Stages):

- 5-class cartilage compartment segmentation (V-Net)
- 11-class segmentation with original image + 5-class input
- ROI extraction with volumetric bounding boxes
- 17 3D CNNs → 52 probabilities for all lesion types


Recon = Reconstruction



# Multi-Plane Efficacy and Grading

## Four-Grade ACL Injury Classification

- **Grading System:**
  - **Grade 0:** No abnormalities
  - **Grade I:** Good continuity, intact contour, <50% damage (Small signal patches/streaks, slightly thickened)
  - **Grade II:** Poor continuity, some fibers visible,  $\geq 50\%$  damage (Local/diffuse thickening, incomplete edges, high signal)
  - **Grade III:** Complete rupture, broken continuity (Displaced/bent ends, clumped ligament, unclear boundaries)
- Model Architecture:



# Category 4: Advanced Architecture Approaches



## Innovative Network Designs for Enhanced Performance

- **Focus:** Developing **high-performance models** with **minimal training data**
- **Key Innovations:**
  - **Attention mechanisms** for capturing **multi-scale relationships**
  - **Graph Convolutional Networks** for **spatial dependencies**
  - **Self-supervised contrastive learning**
  - **Hybrid attention modules** (layer, channel, space scales)
- **Clinical Value:**
  - Improved accuracy through better feature representation
  - Faster inference times (down to **0.5s per scan**)
  - Reduced false positives and false negatives

# SGNET

## Selective Group Attention Network

- **Design Philosophy:** Mimics clinical diagnosis workflow
- **Three Key Modules:**
  - 1. **Dual-Scale Data Augmentation (DDA):** Enriches training data on **spatial & layer scales**
  - 2. **Selective Group Attention (SG):** Captures **relationships across layer, channel, and space**
  - 3. **Fusion Module:** Explores **inter-relationships** among perspectives | Achieves final classification
- **Advantage:** Multi-scale attention mechanism improves detection of subtle ACL tears



The pipeline of data processing.  
step 1. random crop; step 2.  
normalization; and step 3. erasing  
or mixup


# KneeXNet

## Graph Convolutional Networks for Spatial Dependencies

- **Novel Approach:** Represents knee MRI as a graph
  - **Nodes:** Key anatomical landmarks
  - **Edges:** Spatial connections between these landmarks
- **Three Main Components:**
  - **Graph Construction:** Extracts features from local patches
  - **GCN Layers:** Learn hierarchical features
  - **Multi-Scale Fusion:** Combines features across resolutions
- **Key Innovation:**
  - Self-supervised contrastive learning
  - Positive/negative MRI patch pairs
  - Enhances discriminative/robust representations
- **Performance:**
  - Rapid: **0.5s** per MRI scan
  - Ablation studies confirm **GCN layers essential for spatial dependencies**
- **Results:**
  - AUC scores: **0.985**
  - ACL tears detection: **0.972**
  - Meniscal tears: **0.968**

# KneeXNet

## Graph Convolutional Networks for Spatial Dependencies



The overall architecture of KneeXNet

# ACGNN

## Attentive Graph NN with CRF and Meta-Learning

- **Classification Task:** 6 grades (Types I-V based on tear location + Normal)


The tear location indicates the **distal remnant length** as percentage of the ligament length

Tear Type	Type I	Type II	Type III	Type IV	Type V
Location %	> 90	75-90	25-75	10-25	<10

- **Three Key Innovations:**
  - Lightweight Feature Embedding (**knowledge distillation**):
    - Better accuracy-compute trade-off
  - **CRF Integration:**
    - **Refines sample affinities** in GNN layers
  - **Meta-Learning Strategy:**
    - Fast/accurate classification
    - CRF loss addition boosts accuracy >1%
- **Advantage:** Precise tear location classification guides surgical planning

# ACGNN

## Attentive Graph NN with CRF and Meta-Learning



Overview of the pipeline of the proposed DL framework

# Category 5: Data-Efficient & Few-Shot Learning


## Overcoming Data Scarcity Challenges

- **Focus:** Novel architectural innovations including attention mechanisms, graph networks, and specialized feature extraction
- **Key Challenge:**
  - Medical imaging datasets are **expensive and time-consuming to annotate**
  - **Limited samples for rare pathologies**
  - Need for **specialized, site-specific models**
- **Solution Approach:**
  - **Few-Shot Learning (FSL)** frameworks
  - **Metric learning with optimized embedding spaces**
  - **Transfer Learning (TL)** from larger datasets
- **Clinical Value:**
  - Enables deployment in data-scarce clinical settings
  - Faster model development for specialized applications
  - Reduces annotation burden on radiologists

# MedNet-FS

## 3D Few-Shot Learning Framework


- **Goal:** Substantially reduce data dependency for specialized medical imaging tools
- **Technical Approach:**
  - Few-Shot Algorithm integrated with generalized end-to-end (**GE2E**) loss
  - Optimizes embedding space:
    - ↓ **Intra-class distances** (similar cases closer)
    - ↑ **Inter-class distances** (different cases farther)
    - → More **discriminative features**
- **FSL Concept:**
  - Data-efficient approach for recognizing unseen classes
  - **k-shot task:** learns from k examples per class
- **Performance:**
  - **AUC of 0.76** for ACL tear classification, using only **k=40 samples**
- **Impact:** Demonstrates viability of high-performance models with minimal data



(A) traditional Cross-Entropy Loss; (B) the GE2E Loss, designed to optimize the embedding space

# MedNet-FS

## 3D Few-Shot Learning Framework



Supervised training process with MedNet (ResNet101)

# Category 6: Multi-Plane & Multi-Task Models


Comprehensive Analysis Through Integration

- **Focus:** Leveraging multiple MRI planes and simultaneous task execution for comprehensive diagnosis
- **Key Approaches:**
  - Multi-plane fusion (sagittal, coronal, axial)
  - Multi-task learning (segmentation + classification)
  - Semi-supervised learning (pseudo-labeling)
  - Cross-sequence attention mechanisms
- **Clinical Value:**
  - Mimics radiologist workflow (reviewing multiple planes)
  - Simultaneous diagnosis of multiple abnormalities
  - Improved accuracy through complementary information
  - Reduced false negatives by cross-validating findings

# TripleMRNet

## Multi-Plane Feature Fusion


- **Capability:** Analyzes 7 flexible combinations of sagittal, coronal, and axial inputs
- **Design Philosophy:**
  - Mimics clinical workflow of synthesizing multiple planes
  - **Feature-level fusion via concatenation**
- **Efficiency Metrics (Priority):**
  - F1 Score: Balancing precision ( $\downarrow$  false positives)
  - Recall: Avoiding missed injuries ( $\downarrow$  false negatives)
- **Performance Results:**
  - ACL Tear Detection:
    - Three-plane model: Best performance
      - **Accuracy: 0.925 | F1 Score: 0.919**
  - Meniscal Tears:
    - **Axial plane** surpassed **sagittal plane**
- **Key Finding:** Multi-plane integration provides complementary information



# TripleMRNet

## Multi-Plane Feature Fusion

The **orange border** indicates the final convolutional layer computing the **CAM scores** and **generate heatmaps**.



# CoPAS

## Co-Plane Cross-Sequence Attention Network

- **Scope:** Multi-task model for simultaneous diagnosis of **12 knee abnormalities**
- **Technical Innovation:**
  - Generates **synthetic cross-plane volumes**
  - **Integrates via attention mechanism**
  - Captures **spatial correspondence and inter-sequence dependencies**
- **Performance:**
  - Strong generalizability across multiple centers
  - Surpasses SOTA models across all 12 classes
  - Cross-dataset evaluation demonstrates robustness
- **Clinical Impact:**
  - Boosts diagnostic performance for junior AND senior radiologists
  - Serves as effective clinical assistant
  - Class activation maps provide interpretability
- **Advantage:** Comprehensive knee assessment in single forward pass






Illustration of class activationmap  
on a patient with ACL

# CoPAS


## Co-Plane Cross-Sequence Attention Network



The pipeline of CoPAS

# CoPAS


## Co-Plane Cross-Sequence Attention Network



# SDMT

## Spatial Dependence Multi-Task Transformer


- **Tasks:** Joint 3D knee MRI segmentation + landmark localization
- **Architecture Innovations:**
  - Hybrid Multi-Head Attention:
    - **Inter-task attention:** Spatial dependence between tasks
    - **Intra-task attention:** Correlation within single task
  - Shared encoder for resource efficiency
  - **Dynamic weighted multi-mask loss:** Balanced optimization between tasks
- **Advantage:** Tasks mutually enhance each other through spatial dependencies



The knee segmentation objectives and landmark localization objectives.

# SDMT

## Spatial Dependence Multi-Task Transformer



### Three-Stage Pipeline:


- **Shared Encoder** → Multi-scale feature extraction
- **SDMT Module** → Captures spatial dependencies + corrects features between tasks
- **Dual Decoders** → Independent segmentation and landmark outputs

The overall network structure of the proposed method

# DCLU-Net

## Dual-Module System with Semi-Supervised Learning

- **Architecture:** Dual-module U-Net-based system
  - **Module 1:** Segmentation (3D-UNet based)
  - **Module 2:** Classification (Conv-3D + Linear layers)
- **Semi-Supervised Learning:**
  - Teacher model generates pseudo masks
  - Applied to 150 unlabeled samples
  - Reduces manual segmentation needs
- **Novel Feature:** Radiomic Features Integration
  - Extracted from segmented ACL masks
  - Backpropagated to classification module
  - Enhances classification ability
- **Advantage:** Leverages unlabeled data to improve performance




Illustrative examples of the model (with pseudo data) for segmentation of three classes and contrast-enhanced visualization of ACL region.

# DCLU-Net

## Dual-Module System with Semi-Supervised Learning

Overview of semi-supervised segmentation framework




(A) DCLU-Net model with both segmentation & classification modules with extra radiomic features

(B) The segmentation: based on the 3D-UNet  
The classification module combined with two Conv-3D layers and two linear layers

# One-Stop Detection

YOLOv5-Based Integrated System


- **Two Main Components:**
  - **ACL Localization (YOLOv5m)**
  - **ACL Multi-Classification:**
    - **Squeeze-and-Excitation (SE) blocks**
    - Integrates multi-slice data via attention weights
- **Outputs:**
  - ACL localization bounding box
  - Classification result
  - Probability plot for confidence visualization
- **Advantage:** 'One-stop' solution from localization to diagnosis




Sagittal MRI views of (A) a correctly classified knee with a partial tear ACL and (B) its ACL localization with, (C) probability plot

# One-Stop Detection


## YOLOv5-Based Integrated System



(A) ACL localization model



(B) ACL injury classification



Overview of the one-stop detection model pipeline

# Conclusion

- Major Achievements:
  - Segmentation: 1-2 hours → <1 second per case
  - Multi-class grading approaching clinical standards
  - Few-shot learning enabling data-efficient deployment
  - Multi-plane models mimicking radiologist workflow
- Clinical Integration Progress:
  - Improved interreader agreement (Cohen's  $\kappa$ : 0.42 → 0.61)
  - Real-time processing capabilities (0.5s per scan)
  - Multi-task systems diagnosing 12+ abnormalities simultaneously

**Thank you for your time and attention!**